E1

User Manual

breo®

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ENGLISH

User Manual

Thank you for purchasing the breo products.
Before using this product,
please read this manual carefully
for convenient use.

Contents

Product Overview and Components	25
How to Use	26
Product Specifications and Troubleshooting	27
Maintenance	28
Cautions	29
Guidance on Product Replacement/Refund	30
Guidance on Quality Assurance	31
Product Warranty	32

Product Overview and Components

Product Overview

E1 will relax your eyes through smart air-driven acupressure, vibration and heating. Its foldable, luxurious design along with its music playing feature will ease your mind and body. Furthermore, its soft material where your eyes are contacted will give you a best fit.

Product Components

E1, User Manual, Charging Adapter

E1

User Manual

Charging Adapter

How to Use

Preparations before Use

1. Make sure that the battery has been fully charged.
2. You should always charge the battery with a reliable charging adapter that comes with the Product.
3. Wear the Product and turn it on. This step is designed to prevent any air bag damage from too much air pressure.

Charging the Product

1. open the cover of the Product's charging terminal.
2. Plug the charging adapter's DC plug into the charging terminal.
3. The battery level icon on the status indicator blinks while charged.
4. The Product will automatically stop while it is being charged.
5. It takes about 2 to 3 hours to charge the Product.

Power / Control

1. If you press and hold the ON/OFF button for 2 seconds, the relaxer will be turned on.
2. If you long-press the ON/OFF button, it will be turned off.
3. Mode selection: If you short-press the ON/OFF button after turning the relaxer on, mode will change in a preset sequence (SLEEP mode → MED mode → HARD mode).

	SLEEP Mode	MED Mode	HARD Mode
Time	10 minutes	15 minutes	20 minutes
Air Pressure	O	O	O
Heating	O	O	O
Vibration	X	X	O
Music	X	O	O

Product Specifications and Troubleshooting

Product Specifications

Rated Voltage	DC 5V 1A
Power Consumption	5W
Dimensions / Weight	(W)208mmX(L)76mmX(H)106mm / 318g
Operating Environment	Temperature : 5°C ~ 40°C / Humidity : ≤80% Air pressure : 860hpa ~ 1060hpa

Troubleshooting

If you experience any problem using the Product, please check the methods specified below. If the problem persists, turn the Product off and send it to the service center. Never disassemble or repair the Product on your own.

Problems	Possible Causes	Solution
If the Product cannot be turned on	Low or dead battery	Charge the battery.
If air pressure or vibration is too weak	Low or dead battery	Charge the battery.
If the Product cannot be charged	Poor charger connection	Make sure the charger has been properly connected.

Maintenance

Storage precautions

Keep the product in a safe place to avoid exposure to moisture and contamination. If the storage period exceeds 6 months, check the operation of product and turn on the power for about 30 minutes. Avoid the direct sunlight and corrosive gases.

Transporting precautions

This product is designed for seismic resistance and can be transported by plane, train, car, and ship. Avoid snow, rain, and impact as much as possible and put it in the right direction.

Troubleshooting

If a malfunction occurs during operation, check the main battery first. If the trouble cannot be solved, turn off the power and send it to our A/S service center. Do not self-disassemble or repair. If the voltage exceeds the rated range of product, do not turn on the product. Check the cables of power supply periodically.

Cautions

※ Please read the following cautions carefully and accordingly use the Product.

1. Do not use the product after you drink, when you don't feel good, or while you are driving.
2. Any person with a stroke, vascular disease, hypoglycemia or eye disease (e.g. glaucoma, cataract, retinal detachment, etc.) is not allowed to use the Product.
3. Any child aged under 7, any immobile patient, or any physically or mentally challenged person is not allowed to use the Product.
4. Any person with a pacemaker or electronic prosthesis, or any pregnant woman is not allowed to use the Product.
5. Do not wear contact lenses or glasses during the Product.
6. Any person who has had any eye surgery (e.g. refractive surgery, plastic surgery, etc.) is not recommended to use the Product. If such person still intends to use the Product, he or she shall consult with his or her doctor.
7. It is recommended to use the Product with your eyes closed. Once you're done using the Product, do not immediately open your eyes but slowly open them in a few seconds.
8. The Product does not substitute any diagnosis or treatment by a medical specialist. If you want to confirm that you have no health issue for using the Product, please consult with your doctor.
9. If a body part to be treated with the Product has a rash or infection, or has undergone a surgery or treatment, you cannot use the Product.
10. The Product has the heating function. Therefore, any person who is insensitive to heat should take extra care.
11. If you experience any problem or displeasure while using the Product, stop and turn off the Product immediately, and perform relevant inspection.

Cautions

12. You should only use a reliable charging adapter that comes with the Product and take extra care for electrical safety.
13. Wear the Product and turn it on. While you are wearing the Product, you are not supposed to walk or run. Take extra care for electrical safety.
14. The Product is a household product. Therefore, it is strictly prohibited from using the Product for any commercial, treatment or diagnostic purpose.
15. Use a sterilized cotton dipped in neutral detergent when cleaning the Product. Never use water as cleaning liquid.

Guidance on Product Replacement / Refund

1. If a critical repair due to any performance- or function-related defect under a normal operating environment is required within 10 days from the date purchased, you can request for a product replacement or refund with your purchase receipt.
2. If a critical repair due to any performance- or function-related defect under a normal operating environment is required within 1 month from the date purchased, you can request for a product replacement or repair. However, such request can only be valid if the Product has been damaged under a normal operating environment. In addition, the Product can be replaced, refunded or repaired only if the user manual and all other components have maintained their original conditions before unpacked.
If you intend to make such request above or have any related inquiry, please contact your buyer or our customer center (1588-2336). Batteries are not applicable to any warranty service.

Guidance on Quality Assurance

Warranty information

With the product warranty, you can easily receive various maintenance services. Only one warranty card is issued for each product. Therefore, you should take extra care not to lose the card since it can be really cumbersome to reissue a new card.

※The warranty period is one year from the date of purchase.

Warranty conditions

1. With the product warranty, you can receive various repair services for one year from the date purchased. However, such services may not be properly provided if the serial number on the Product is not legible, or if such serial number differs from the one registered on the product warranty. (Batteries are not applicable to such services.)
2. For any defect due to a user's fault (e.g. use of wrong power, accidental failure, etc.) or a repair service by a non-authorized vendor, no warranty service can be provided.
3. Additional expense may incur for any supplementary services not specified in the product warranty (e.g. delivery service, etc.)
4. For any defect due to a user's fault or negligence, a disassembly service by an unqualified expert, no product warranty, or a natural disaster, no warranty service can be provided.
5. When requesting for any quality assurance service, please enclose your product warranty. We are responsible for collecting, treating and recycling any returned or disposed product.
6. If the Product has been used for any commercial purpose (e.g. business use, public utilities, etc.) other than household uses, no warranty service can be provided.
7. For any defect due to a natural disaster (e.g. typhoon, flood, earthquake, etc.), no warranty service can be provided.

※ Note: Since quality assurance services get continuously upgraded, detailed regulations may be subject to change without any customer notification.

Product Warranty

Name of Product : _____

Serial Number : _____

Date Purchased : _____

Place Purchased : _____

Name of Customer : _____

Contact : _____

Address : _____

Guidance on Product Replacement/Refund

1. If a critical repair due to any performance- or function-related defect under a normal operating environment is required within 10 days from the date purchased, you can request for a product replacement or refund with your purchase receipt.
2. If a critical repair due to any performance- or function-related defect under a normal operating environment is required within 1 month from the date purchased, you can request for a product replacement or repair. However, such request can only be valid if the Product has been damaged under a normal operating environment. In addition, the Product can be replaced, refunded or repaired only if the user manual and all other components have maintained their original conditions before unpacked. If you intend to make such request or have any related inquiry, please contact your buyer or our customer center (1588-2336).

Product Warranty

Warranty information

With the product warranty, you can easily receive various maintenance services. Only one warranty card is issued for each product. Therefore, you should take extra care not to lose the card since it can be really cumbersome to reissue a new card.

※ The warranty period is one year from the date of purchase.

Warranty conditions

1. With the product warranty, you can receive various repair services for one year from the date purchased. However, such services may not be properly provided if the serial number on the Product is not legible, or if such serial number differs from the one registered on the product warranty. (Batteries are not applicable to such services.)
 2. For any defect due to a user's fault (e.g. use of wrong power, accidental failure, etc.) or a repair service by a non-authorized vendor, no warranty service can be provided.
 3. Additional expense may incur for any supplementary services not specified in the product warranty (e.g. delivery service, etc.)
 4. For any defect due to a user's fault or negligence, a disassembly service by an unqualified expert, no product warranty, or a natural disaster, no warranty service can be provided.
 5. When requesting for any quality assurance service, please enclose your product warranty. We are responsible for collecting, treating and recycling any returned or disposed product.
 6. If the Product has been used for any commercial purpose (e.g. business use, public utilities, etc.) other than household uses, no warranty service can be provided.
 7. For any defect due to a natural disaster (e.g. typhoon, flood, earthquake, etc.), no warranty service can be provided.
- ※ Note: Since quality assurance services get continuously upgraded, detailed regulations may be subject to change without any customer notification.