# NP 1 air

User Manual

breo®

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User Manual

breo®

Thank you for purchasing the breo products.  
Before using this product,  
please read this manual carefully  
for convenient use.

## Contents

---

Product Overview and Components	30
How to use	31
Specifications and Troubleshooting	33
Maintenance	34
Precautions	35
Replacement and refund	36
Warranty Guide	37
Product Warranty	38

## Product Overview and Components

---

### Product Overview

NP1 air massages your neck closely wrapping your neck snugly with soft, light and fully filled air.

### Product Components

NP1 air, Strap, Hard case, User Manual, Charging adaptor

NP1 air

Strap

Hard case

User Manual

Charging adaptor

## How to use

---

### Preparation before use

1. Please check the product before use.
2. Make sure you have all of the components.
3. Make sure that the charging cable is not damaged.

### Charging

1. Open the charging cover on the side of the main unit and charge it using the charging cable supplied.
2. Red light blinks when battery is low. A red LED is lit during charging and a green LED is lit when charging has been completed.
3. The product will not work during charging. If you attempt to charge the battery while operating, the self-protection program will be executed to stop operation automatically.
4. It can be used for 5 hours based on full charge.

### Wearing

1. Attach the supplied strap to the product.
2. Use the product as if wearing a neck pillow.

## How to use

---

### How to use

1. Make sure the product is fully charged.
2. Make sure the air tube is completely closed.
3. Press and hold the power button to turn on the power and the air will be automatically filled when the product runs out of air. (About 110 seconds)
4. When the air is completely filled, the product will start operation automatically.
5. The product usage time is automatically set to 15 minutes after starting massage.
6. You can select one of 2 rotation directions according to your preference.
7. You can have stronger massage by connecting and pulling the strap.
8. Press and hold the power button to stop the operation and you can use it like a regular neck pillow.
9. If you feel that the product is running out of air when using your neck pillow, you can fill it after turning on the product again.
10. Open the air tube after using the product to remove the air completely, and the volume will decrease. If necessary, store it in a convenient hard case.

### Mode Selection

Mode 1	Mode 2
Right direction (Default)	Reverse direction

## Specification and Troubleshooting

---

### Specifications

Rating Voltage	DC 5V 1A		
Battery	2300mAh / 3.7V / Lithium polymer	Power consumption	5W
Materials	ABS / Sponge / Cloth	Charging Terminal	C-Type
Size/Weight	(W)310mm X (L)310mm X (H)120mm / 650g		

### Troubleshooting

If you have any problem in using the product, please check the instructions below. If this does not work, turn off the power and send it to the service center. Please note that if you are not an expert, do not disassemble or repair the product personally since we may not provide a further service.

Problem	Possible Cause	Solution
Power is not turned on	Battery is low or discharged	Charge the battery
Rotation gets weak	Battery is low or discharged	Charge the battery

## Maintenance

---

### Storage precautions

Keep the product in a safe place to avoid exposure to moisture and contamination. If the storage period exceeds 6 months, check the operation of product and turn on the power for about 30 minutes. Avoid the direct sunlight and corrosive gases.

### Transporting precautions

This product is designed for seismic resistance and can be transported by plane, train, car, and ship. Avoid snow, rain, and impact as much as possible and put it in the right direction.

### Troubleshooting

If a malfunction occurs during operation, check the main battery first. If the trouble cannot be solved, turn off the power and send it to our A/S service center. Do not self-disassemble or repair. If the voltage exceeds the rated range of product, do not turn on the product. Check the cables of power supply periodically.

## Precautions

---

※ Please read the following cautions carefully and accordingly use the Product.

1. Do not use the product if you have any of the following symptoms (surgery area, skin irritation, erythema, inflammation, heart disease, acute disease, skin disease, malignant tumor, blood disease, bleeding, osteoporosis, soft tissue disease, stroke, chronic illness, pregnant woman or woman in menses) and use the product after consulting your doctor.
2. If you have a physical or mental disability and are athletic or mentally retarded, please use the product after consulting your doctor.
3. This product is not a medical device, so please do not use the product for the purpose of diagnosis and treatment.
4. Do not use the product in hunger or after overeating, drinking or excessive exercise.
5. Do not use the product while driving.
6. Please use the charger correctly when charging and be careful for electric safety.
7. Avoid environments such as water or fire which may damage the product.
8. Do not drop or damage the product with a sharp object.
9. The minor shall use the product together with a guardian.
10. If any problem occurs during use, stop operation immediately and contact Customer Center (Do not disassemble the product arbitrarily).
11. If you feel uncomfortable during use, please stop operation immediately.

## Guidance on Product Replacement / Refund

---

1. If a critical repair due to any performance- or function-related under a normal operating environment is required defect within 10 days from the date purchased, you can request for a product replacement or refund with you purchase receipt.
2. If a critical repair due to any performance- or function-related defect under a normal operating environment is required within 1 month from the date purchased, you can request for a product replacement or repair. However, such request can only be valid if the Product has been damaged under a normal operating environment. In addition, the Product can be replaced, refunded or repaired only if the user manual and all other components have main tained their original conditions before unpacked.

If you intend to make such request above or have any related inquiry, please contact your buyer or our customer center (1588-2336). Batteries are not applicable to any warranty service.

## Guidance on Quality Assurance

---

### Warranty information

With the product warranty, you can easily receive various maintenance services. Only one warranty card is issued for each product. Therefore, you should take extra care not to lose the card since it can be really cumbersome to reissue a new card. ※ The warranty service is valid for one year from the date of purchase.

### Quality Assurance Conditions

1. With the product warranty, you can receive various repair services for one year from the date purchased. However, such services may not be properly provided if the serial number on the Product is not legible, or if such serial number differs from the one registered on the product warranty.  
(Batteries are not applicable to such services.)
2. For any defect due to a user's fault(e.g. use of wrong power, accidental failure, etc.) or a repair service by a non-authorized vendor, no warranty service can be provided.
3. Additional expense may incur for any supplementary services not specified in the product warranty.(e.g. delivery service, etc.)
4. For any defect due to a user's fault or negligence, a disassembly service by an unqualified expert, no product warranty, or a natural disaster, no warranty service can be provided.
5. When requesting for any quality assurance service, please enclose your product warranty. We are responsible for collecting, treating and recycling any returned or disposed product.
6. If the Product has been used for any commercial purpose(e.g. business use, public utilities, etc.) other than household uses, no warranty service can be provided.
7. For any defect due to a natural disaster(e.g. typhoon, earthquake, flood, etc.), no warranty service can be provided.

Note: Since quality assurance services get continuously upgraded, detailed regulations may be subject to change without any customer notification.

## Product Warranty

---

Name of Product : \_\_\_\_\_

Serial Number : \_\_\_\_\_

Date Purchased : \_\_\_\_\_

Place Purchased : \_\_\_\_\_

Name of Customer : \_\_\_\_\_

Contact : \_\_\_\_\_

Address : \_\_\_\_\_

### Guidance on Product Replacement/Refund

1. If a critical repair due to any performance- or function-related under a normal operating environment is required defect within 10 days from the date purchased, you can request for a product replacement or refund with your purchase receipt.
2. If a critical repair due to any performance- or function-related defect under a normal operating environment is required within 1 month from the date purchased, you can request for a product replacement or repair. However, such request can only be valid if the Product has been damaged under a normal operating environment. In addition, the Product can be replaced, refunded or repaired only if the user manual and all other components have main tained their original conditions before unpacked If you intend to make such request above or have any related inquiry, please contact your buyer or our customer center(1588-2336). Batteries are not applicable to any warranty service.

# Product Warranty

---

## Warranty information

With the product warranty, you can easily receive various maintenance services. Only one warranty card is issued for each product. Therefore, you should take extra care not to lose the card since it can be really cumbersome to reissue a new card.

※ The warranty service is valid for one year from the date of purchase.

## Quality Assurance Conditions

1. With the product warranty, you can receive various repair services for one year from the date purchased. However, such services may not be properly provided if the serial number on the Product is not legible, or if such serial number differs from the one registered on the product warranty.  
(Batteries are not applicable to such services.)
2. For any defect due to a user's fault(e.g. use of wrong power, accidental failure, etc.) or a repair service by a non-authorized vendor, no warranty service can be provided.
3. Additional expense may incur for any supplementary services not specified in the product warranty.(e.g. delivery service, etc.)
4. For any defect due to a user's fault or negligence, a disassembly service by an unqualified expert, no product warranty, or a natural disaster, no warranty service can be provided.
5. When requesting for any quality assurance service, please enclose your product warranty. We are responsible for collecting, treating and recycling any returned or disposed product.
6. If the Product has been used for any commercial purpose(e.g. business use, public utilities, etc.) other than household uses, no warranty service can be provided.
7. For any defect due to a natural disaster(e.g. typhoon, earthquake, flood, etc.), no warranty service can be provided.

Note: Since quality assurance services get continuously upgraded, detailed regulations may be subject to change without any customer notification.